




DR. Chip INTRODUCTION

**Bringing Innovative Products and Total Solutions to
the Modern Molecular Diagnostics**



What is DR. ?

DR. =DNA+RNA

=Doctor


=Expert

DR. Chip aims to use DNA and RNA
to give a new tool for diagnosis

The Needs of a Biochip


Biochip vs Multiplex PCR

Numerous reasons including following,




Good multiplex PCR and electrophoresis but difficult to identify correct fragments

PCR




- *Streptococcus* spp.
- *Strep. agalactiae*
- *Strep. bovis*
- *Strep. dysgalactiae*
- *Strep. uberis*
- *E. coli*
- *Staphyl. aureus*
- Positive control
- Negative control
- Blank



Easy-recognized pattern with sequence validation by a chip

Design Flow Chart

- Search target genes (Gene bank)
- Primer design, testing, and validation
- Probe design, testing, and validation



Kit contains

- Extraction Reagents*
- Primer Mixture
- DNA Polymerase
- DR. Hyb Buffer
- Blocking Reagent
- Wash Buffer
- Strep-AP
- NBT/BCIP
- Detection Buffer





Products

FOOD Related Products

- DR. Food Chip Series – Foodborne pathogens
- DR. Milk Chip – Mastitis pathogens
- DR. WIT Chip – Marine Shrimp Viral Diseases

Products

Human Related Products

- DR. RV Chip - Respiratory viruses, including SARS
- DR. EV Chip - Enterovirus
- DR. HPV Chip - Human Papillomavirus
- DR. RB Chip - Respiratory bacteria
- DR. TB Chip - Mycobacterium tuberculosis
- DR. HBV Chip - Hepatitis Type B & Lamivudine drug resistant mutations.

Coming up...


DR. Livestock DNA Chip

DR. Orchid virus Chip

DR. Porcine Chip


platform

- Polychip
- New 3 well
- C8 strips
- 96 well



Platform Characteristic

- Low noise
- Easy to use
- Cost saving



DB Chip polymer chip (low noise)





DR. AiM™ Technology Arrays-in-Microplate

High Throughput Screening Array
Array Reader and
Automatic Analysis System

DR. AiM™


Microplate-based Diagnostic Array

◆ DR. AiM™ - standard 96-well plate, which is pre-spotted with specific probes at the bottom of the well to capture gene(s) of interest.



Key Components

- Microarray Plate
- Mini Oven
- Fluidic Station
- Biochip Reader
- Biochip Analysis Software



Semi-automatic system

**DNA
Extraction**



PCR




**Hybridization
oven**



reading



**Color
developing**



DR. AiM™ Reader

- ◆ **DR. AiM™ Reader**- a sophisticate device that capture image of DR. AiM™



Designed to accommodate

1. standard 96-well microtiter plate
2. DR. PolyChip (3 wells)

DR. PolyChip (slide format, 3 well / slide)



96-well microtiter



DR. AiMsoft™

DR. AiM Reader

DR. AiM Reader

SSS Library

File Setup

HW Setup

Print

Data Analysis

Admin Login

Help

Exit

DR. AiM Reader

Processing finished!

DR. AiM Reader

Show / Modify Result

Probe Diagram

Probe Name
Hybridization Positive control
Hybridization Negative Control
PCR Positive Control
Staphylococcus aureus
Escherichia coli
Salmonella spp.

Feature Diagram for Manual Comparison

= Negative Controls

Feature	Selected ?
Staphylococcus aureus	
Escherichia coli	
Salmonella spp.	Yes

Result

Well Position : E7

Result : Positive

You may modify the result below :

Positive Negative Invalid

Feature	Present ?
Staphylococcus aureus	
Escherichia coli	
Salmonella spp.	Positive

Cancel Apply

Delete Image


User Mode : Regular User

Total # of Wells = 96 # of Positive = 72 # of Negative = 24 # of Invalid = 0

Positive Manual Invalid

BioChip Analysis Software

- User-friendly programmable interface
- Imaging and Analysis
 - Background noise reduction
 - Template Library Directory Setup
 - Data Storage and Output in Excel
 - Default Setup



DR. Fluidic soft™

DR.FluidicSoft - [Food940121]

File View Instrument Options Window Help

Procedures Cycles Plates Sequence

Volume Settings for Prime & Rinse

Rinse : 30 ml Wash1 : 15 ml Wash2 : 15 ml Wash3 : 15 ml

Sequence of Procedures

Start at Index : 1 Number of Strips : 0 (All Strips)

Index	Procedure / Prime&Rinse	Time Gap (HH : MM : SS)
1	Rinse	00 : 00 : 00
2	Wash1	00 : 00 : 00
3	5 WASH	00 : 00 : 00
4	Rinse	00 : 00 : 00
5	Wash2	00 : 00 : 00
6	BLOCK	00 : 00 : 00
7	Rinse	00 : 30 : 00
8	Wash1	00 : 00 : 00

Add Step Insert Step Modify Step Delete Step Delete All Execute

Ready User Mode : Regular User

開始 DR.FluidicSoft - [Foo... CH 下午 04:01



Fluidic Station

- A semi-automatic fluidic handling system
 - Eliminate human error
- User-friendly programmable interface
- For 96-well and 8 strips microplate
- Compact in size
 - 6.2 kg
 - 28*22*42 (cm)



Mini Oven

- Accurate Temperature Control
 - Range: Room Temp. ~ 85°C
- Vibration to enhance Hybridization
 - 60 Hz
- Compact in Size
 - 6.2 kg
 - 26*30*17.5 (cm)
- Fast Heat Up





DR. Food™ Kit series

Identify maximum 10 major foodborne pathogens in one test

Staphylococcus aureus, Salmonella spp. E. coli
Yersinia enterocolitica, Bacillus cereus,
Clostridium perfringens, Listeria monocytogenes
Shigella spp., Vibrio spp., Campylobacter



◆ Application-

- **Detection of common foodborne bacteria**
- **Assist with quality control of food processing**
- **Rapid method for foodborne pathogen disease**

◆ Process

Food samples →

enrichment overnight →

extract DNA → DNA amplification →

hybridization → color indication



Kit name	Detecting Bacteria	note
DR. Food – 9 Kit	Staphylococcus aureus, Salmonella spp. E. coli Yersinia enterocolitica, Bacillus cereus, Clostridium perfringens, Listeria monocytogenes Shigella spp. Vibrio spp.	Universal pathogens detection
DR. Food –10 Kit	Staphylococcus aureus, Salmonella spp. E. coli Yersinia enterocolitica, Bacillus cereus, Clostridium perfringens, Listeria monocytogenes Shigella spp. Vibrio spp. Campylobacter	Universal pathogens detection
DR. Food – S Kit	Staphylococcus aureus, Salmonella spp. E. coli, Listeria monocytogenes Vibrio cholerae.	QC for Sea food material process



DR. Food™ Kit

● Designing and Probe Allocation



●	Hybridization Positive Control
●	<i>Staphylococcus aureus</i>
●	<i>Salmonella</i> spp.
●	<i>Escherichia coli</i> - <i>Shigella</i> spp.
●	<i>Listeria monocytogenes</i>
●	<i>Bacillus cereus</i>
●	<i>Yersinia enterocolitica</i>
●	<i>Vibrio</i> spp.
●	<i>Clostridium perfringens</i>
●	PCR positive control
●	Negative control



Performance

Detection Limit (cfu/mL)

Pathogen	Sensitivity	Pathogen	Sensitivity
<i>Staphylococcus aureus</i>	$10^9 \sim 10^4$	<i>Listeria monocytogenes</i>	$10^8 \sim 10^4$
<i>Escherichia coli</i>	$10^8 \sim 10^4$	<i>Salmonella</i> spp.	$10^9 \sim 10^4$
<i>Yersinia enterocolitica</i>	$10^8 \sim 10^4$	<i>Shigella</i> spp.	$10^8 \sim 10^4$
<i>Bacillus cereus</i>	$10^8 \sim 10^3$	<i>Vibrio</i> spp.	$10^9 \sim 10^4$
<i>Clostridium perfringens</i>	$10^8 \sim 10^3$		







Mastitis Factors and Diagnosis

- The most costive cow disease
- According to US National Mastitis Council
 - US\$ 184 per cow
 - US\$ 184,400 per 100 cow herd
 - US\$ 1.7 billion in US
- To detect and identify the pathogens
 - Choosing the right treatment
 - Increasing milk production and milk quality





Sample preparation



- Hyb Positive Control
- E. coli
- *Streptococcus spp.*
- *Staphyl. aureus*
- *Strep. agalactiae*
- *Mycoplasma bovis*
- *Strep. bovis*
- PCR Positive Control
- *Strep. dysgalactiae*
- Blank
- *Strep. uberis*

remove fat and
nucleic acid purification






DR. Milk Chip Diagnostic Process



Colorimetric development

Nucleic acid amplification

Hybridization




Bacteriology Culture vs. DR. Milk Chip

Bacteria species	Bacteriology culture	DR. Milk Chip
<i>Streptococcus</i> spp.	At least 2 days	Detect 7 mastitis-causing bacteria simultaneously in 6 hours
<i>Strep. uberis</i>	5 days	
<i>Strep. bovis</i>		
<i>Strep. agalactiae</i>		
<i>Strep. dysgalactiae</i>		
<i>Staphyl. aureus</i>	5 days	
<i>E. coli</i>	4 ~ 7 days	
<i>Mycoplasma bovis</i>	7~ 10 days	



Advantages for using DR. Milk Chip :

Parallel Diagnosis :

Detect 7 mastitis-causing bacteria (3 contagious and 4 environmental bacteria) simultaneously

- **Early Diagnosis, Early Treatment :**

Milk separately, decrease the dairy-processing fee and the cost for treatments, as well as habitual abuse of antibiotics by diagnosing in the early stage.


Minimize the loss by segregation, observation and treatments from detecting in the early stage.

- **Other Advantages :**

Increase in milk production and quality, as well as marketing competitiveness by inspecting periodically



DR. Livestock DNA Chip



- Porcine-specific probe
- Bovine-specific probe
- Goat/sheep-specific probe
- Livestock conserve probe
- PCR control
- Hybridization control

Application: DR. Livestock™ Kit is a DNA-based system designed for speciation of pork, beef and goat/sheep meat or detection of MBM in feedstuffs

Sensitivity : The detection limit is 0.1% W/W

Process : The entire operation processes require 6 hours

porcine	Bovine	Goat



Porcine Chip


Purpose: To detect most common virus in Swine.

Target :

- CSFV : Classical Swine Fever Virus
- SVDV : Swine Vesicular Disease Virus
- BVDV : Bovine Viral Diarrhea Virus
- PCV : Porcine Circovirus

Application : Farm swine health monitoring, infection control.

- Cooperation with Council of Agriculture




Future Biodiagnostic Chip

Portable Meter(PDA)




Disposable Chip




DR. RV™ IVD Kit


	Hybridization Positive Control		Parainfluenza III
	Hybridization Negative Control		RSV
	PCR Positive Control		Influenza A
	Coronavirus (SARS)		Influenza B
	Parainfluenza I		Adenovirus
	Parainfluenza II		




INF A




INF B




RSV




Coronavirus (SARS)




PIV 1



PIV 2



PIV 3



Adenovirus



DR. RB™ IVD Kit

Mycoplasma pneumoniae

Chlamydia pneumoniae

Legionella pneumophila



DR. HPV™ IVD Kit

- Detail clarify the quantity and types of HPV when testing with Pap test
- Keep tracking the growing process of cervical cancer definably
- Detecte many different types of HPV accurately !
- Including many type →
HPV high-risk group: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, & 68
HPV low-risk group: 6, 11, 42, 43, & 44



HBV Diagnostic Chip

● Purpose-

- Detect and identify Hepatitis Type B viruses and genotyping A- F.
- Detect Lamivudine drug resistant mutations.

● Process-

- Serum → virus DNA → PCR → hybridization → color indication

● Application-

- Long-term monitoring of HBV treatment and drug therapy

● Collaborating Institute-

- National Taiwan University Medical Center



DR. TB™ IVD Kit

Mycobacterium tuberculosis complex (MTBC)

- *Mycobacterium tuberculosis*
- *Mycobacterium bovis*
- *Mycobacterium bovis BCG*
- *Mycobacterium africanum*
- *Mycobacterium microti*

MOTT or NTM


Mycobacterium other than tuberculosis

Non-tuberculosis *mycobacterium*

Mycobacterium tuberculosis complex screening/drug resistant chip

- **Fast** (Diagnose MTBC infection in 6 hours)
- **High sensitivity** (Lower than 100 bacteria in specimen can be detected)
- **High specificity** (Rare false positive result would be detected with specific probes)
- **High throughput** (Detected in 96 well plate format semi-automatically)
- **Cheap**

Procedure of MTBC screening/drug resistant chip






DR.EV™ IVD Kit

Features:


Size: 1 cm X 1 cm

Probes:


1. Positive & negative controls
2. EV common sequences
3. EV subtype sequences




- Hybridization positive control
- PCR positive control
- EV common sequence
- EV 71
- Cox A16



Enterovirus



EV 71



Cox A16




DR. HBV Chip Kit




Introduction:

1. Detect and identify Hepatitis Type B viruses and genotyping A- F
2. Detect Lamivudine drug resistant mutations (YVDD, YIDD) during HBV therapy.
3. Sensitivity: 10^5 copies/ml
4. Long-term monitoring of HBV treatment and drug therapy
5. Collaborating with National Taiwan University Medical Center


Design of primers



F4 / R58 →
 PCR control →
 15-F / PiKaR →












DR. HBV™ Chip



- Hybridization Positive Control
- Hybridization Negative Control
- PCR Positive Control
- HBV genotype A
- HBV genotype B
- HBV genotype C
- HBV genotype E
- HBV genotype F
- YVDD strain
- HBV common
- YIDD strain


Pattern image examples

				
HBV genotype A	HBV genotype B	HBV genotype C	HBV genotype D	HBV genotype E
				
HBV genotype F	HBV genotype B YVDD	HBV genotype B YIDD	HBV Negative	

DR. HBV™ Chip


HBV type A (YIDD mutant)

10^5 10^6 10^7 copies / ml




HBV type B (YVDD mutant)

10^5 10^6 10^7 copies / ml




HBV type C (non-mutated)

10^5 10^6 10^7 copies / ml



HBV type D (non-mutated)

10^5 10^6 10^7 copies / ml





Assay protocol:

HBV DNA extraction from patient's serum (200 μ l)

QIAamp Blood Mini kit

DNA amplification (PCR method)

Hybridization with chip

Detection with DR.AiM reader



RFLP vs DR. HBV chip

No.	RFLP	DR HBV chip	No.	RFLP	DR HBV chip
R1	B	B	R31	B	B
R2	B	B	R32	B	B
R3	C	C	R33	B	B
R4	B	B	R34	B	B
R5	C	C	R35	C	C
R6	C	C	R36	B	B
R7	B	B	R37	C	C
R8	B	B	R38	B	B
R9	B	B	R39	B	B
R10	C	C	R40	C	C
R11	B	B	R41	B	B
R12	B	B	R42	C	C
R13	B	C	R43	B	B
R14	B	B	R44	B	B
R15	B	B	R45	B	B
R16	B	B	R46	B	B
R17	B	B	R47	B	B
R18	B	B	R48	B	B
R19	B	B	R49	B	B
R20	B	B	R50	B	B
R21	C	C	R51	B	B
R22	B	B	R52	C	C
R23	C	C	R53	C	C
R24	B	B	R54	C	B
R25	C	C	R55	B	B
R26	B	B	R56	C	C
R27	B	B	R57	B	B
R28	C	C	R58	C	C
R29	B	B	R59	B	B
R30	C	C	R60	B	B

RFLP method is from "FEBS Letter" 1999", volume 450, Issue1-2, P66-71. "HBV genotype assignment using RFLP patterns"

DNA Sequencing vs DR. HBV Chip

No.	DNA Sequencing	DR. HBV chip	No.	DNA Sequencing	DR. HBV chip
D1	YIDD	YIDD	D22	YIDD	YIDD
D2	YVDD	YVDD	D23	YVDD	YVDD
D3	YIDD	YIDD	D24	YIDD	YIDD
D4	YVDD	YVDD	D25	YVDD	YVDD
D5	YVDD	YVDD	D26	YIDD	YIDD
D6	YMDD	YIDD	D27	YMDD	YMDD
D7	YMDD	YIDD	D28	YVDD	YVDD
D8	YIDD	YIDD	D29	YIDD	YMDD
D9	YMDD	YMDD	D30	YMDD	YMDD
D10	YMDD	YMDD	D31	YIDD	YIDD
D11	YMDD	YMDD	D32	YMDD	YMDD
D12	YMDD	YMDD	D33	YMDD	YMDD
D13	YMDD	YMDD	D34	YIDD	YMDD
D14	YVDD	YVDD	D35	YIDD	YMDD
D15	YVDD	YVDD	D36	YMDD	YMDD
D16	YIDD	YIDD	D37	YVDD	YVDD
D17	YVDD	YVDD	D38	YVDD	YVDD
D18	YVDD	YMDD		TOTAL	84%(32/38)
D19	YVDD	YVDD		YMDD	83% (10/12)
D20	YIDD	YIDD		YVDD	92%(12/13)
D21	YIDD	YIDD		YIDD	77%(10/13)

HBV⁺ Serum Testing

No	HBeAg	DR HBV chip	Q-PCR (copies/mL)	No	HBeAg	DR HBV chip	Q-PCR (copies/mL)
01	– (0.012)	–	N/A	24	– (0.034)	–	N/A
02	+ (2.697)	+ (B, YMDD)	7,185,000	25	– (0.022)	–	N/A
03	– (0.012)	–	N/A	26	– (0.011)	–	N/A
04	– (0.012)	–	N/A	27	– (0.016)	–	N/A
05	+ (2.598)	+ (B, YMDD)	158,100,000	28	+ (0.196)	+ (B, YMDD)	123,300
06	– (0.022)	–	N/A	29	– (0.015)	–	N/A
07	+ (2.952)	+ (B, YMDD)	302,300,000	30	– (0.045)	–	N/A
08	+ (2.869)	+ (B, YMDD)	19,350,000	31	+ (2.914)	+ (B, YMDD)	3,801,000
09	– (0.013)	–	N/A	32	– (0.020)	–	N/A
10	+ (2.620)	+ (B, YMDD)	232,600,000	33	– (0.018)	–	N/A
11	– (0.014)	–	N/A	34	– (0.013)	–	N/A
12	+ (0.302)	–	N/A	35	+ (2.572)	+ (C, YMDD)	400,700,000
13	– (0.013)	–	N/A	36	– (0.022)	–	N/A
14	– (0.014)	–	N/A	37	– (0.018)	–	N/A
15	– (0.033)	–	N/A	38	– (0.019)	–	N/A
16	+ (2.965)	+ (B, YMDD)	5,071,000	39	– (0.013)	–	N/A
17	– (0.011)	–	N/A	40	– (0.018)	–	N/A
18	– (0.013)	–	N/A	41	– (0.030)	–	N/A
19	– (0.012)	–	N/A	42	+ (2.981)	+ (B, YMDD)	124,600,000
20	– (0.008)	–	N/A	43	– (0.046)	–	N/A
21	– (0.032)	–	N/A	44	– (0.029)	–	N/A
22	– (0.014)	–	N/A	45	+ (3.096)	+ (C, YMDD)	142,000,000
23	– (0.014)	–	N/A				

DR. HPV™ Chip Kit

- Detail clarify the quantity and types of HPV when testing with Pap test
- Keep tracking the growing process of cervical cancer definably
- Detect different types of HPV accurately !




CERVICAL CANCER and HPV

- Widely accepted that human papillomavirus, a sexually transmitted virus, is the major precipitant in the evolution from squamous metaplasia to dysplasia in the transition zone of the cervix
- HPV is the necessary but insufficient cause of cervical cancer

- HPV genomes are found in 99% of premalignant and malignant squamous cell cancer cases worldwide (Walboomers, 2001)
- Presence of high risk HPV (HR-HPV) types (16,18,31,33,34,35,39,45,51,52,56,58,59,66,68,70) are generally recognized as the major risk factor for development of cervical cancer



HPV genotyping



Detect 22 HPV subtypes

High risk :


16,18,31,33,35,39,45,51,52,56,58,59,61,66,70,82

Low risk :

6,11,34,53,54,73

Sensitivity


~10³ copies based on
HeLa cell testing





Testing of HPV chip (I)

- Collaborating with KOO FOUNDATION SUN YAT-SEN CANCER CENTER
- Specimens from Low SIL or ASCUS patients



			ASCUS/L-SIL
Dr.HPV CHIP	Positive	well typing	51
		co-infection	6
	Negative		62
Positive rate%			48

Prevalence

In cervical lesions


- Squamous carcinoma: 95% association
- HSIL/CIN II, III: 75 -95%
- LSIL/CINI: 60%
- ASCUS: 30%
- Adenocarcinoma: 60% association

(ref: Cuzick et al, 1992; Schiffman et al, 1993; IARC, 1995; Olsen et al, 1995)

Testing of HPV chip (II)

- Collaborating with a local reference laboratory center
- Random collection of swab samples

Distribution of HPV subtype



Dr. HPV CHIP	Positive	Well typing	96
		Co-infection	36
	Negative		157
Positive rate (%)			45.7



HPV high/low risk chip

H/L system: 6, 11, 16, 18, 31, 33, 35, 39, 45, 51, 52, 54, 56, 58, 59, 61, 66, 70, 82

