



*Life line lab*

## **PRODUCT SPECIFICATIONS**

PRODUCT:	Epstein Barr virus VCA ANTIGEN
CODE NUMBER:	SL0160
MANUFACTURING:	Manufactured by Lifelinelab s.r.l. Attached General Procedures (Enclose A)
SOURCE:	TPA induced P3HR-1 cell line
FORM:	Frozen Liquid
PURIFICATION:	Virus and viral components suspended in glycine buffer. Cultured virus is sonicated and subject to high speed centrifugation in glycine buffer.
INACTIVATION:	Each batch is inactivated by ultrasonication Note: Although inactivated, please handle with appropriate precautions.
APPLICATIONS:	Final working dilutions must be determined by end user.
STORAGE TEMPERATURE:	Frozen -70°C (avoid repeated thaws)
SHIPPING TEMPERATURE:	Frozen





*Life line lab*

**ENCLOSE A**

**EBV VCA ANTIGEN  
PRODUCTION FLOW SHEET**

P3HR-1 cell line TPA induction



Cell pellet collection



Glycine extraction



Sonication in Glycine Buffer



High speed centrifugation / Supernatant collection

