



PRODUCT SPECIFICATIONS

PRODUCT:	RUBELLA ANTIGEN
CODE NUMBER:	SL0110
MANUFACTURING:	Manufactured by LifeLineLab s.r.l Attached General Procedures (Enclose A)
SOURCE:	Rubella strain - "Therien" growth on B-Vero cells
FORM:	Frozen Liquid
PURIFICATION:	Ultracentrifugation on sucrose gradient.
INACTIVATION:	UV light inactivation.
APPLICATIONS:	Final working dilutions must be determined by end user.
STORAGE TEMPERATURE:	Frozen -70°C (avoid repeated thaws)
SHIPPING TEMPERATURE:	Frozen





Life line lab

ENCLOSE A

**RUBELLA ANTIGEN
PRODUCTION FLOW SHEET**

Rubella strain - "Therien" growth on B-Vero cell line



Collection of supernatants



Purification by ultracentrifugation on sucrose gradient



Isolation of antigen band



Resuspend in carbonate buffer

