



PRODUCT SPECIFICATIONS

- PRODUCT:** Streptavidin coated microplate
- CODE:** MP0100
- SIZE:** 12 strips x 8 wells
- FORM:** Protein adsorbed to plastic and desiccated under vacuum
- PACKAGE:** Sealed under vacuum in aluminium bags containing silica gel.
- APPLICATIONS:** Universal tool for binding biotinilated molecules (e.g. proteins, peptides, polysaccharides, oligonucleotides, DNA fragments).
- STORAGE:** Unopened: Store between +2°C and +8°C
Opened: Put back inside aluminium bag with silica gel, seal with tape and store between +2°C and +8°C
- NOTE:** Not for use in humans or clinical diagnosis. This product is intended for research or manufacturing use only

