



Life line lab

PRODUCT SPECIFICATIONS

PRODUCT:	TOXOPLASMA GONDII ANTIGEN
CODE NUMBER:	SL0120
MANUFACTURING:	Manufactured by Lifelinelab s.r.l. Attached General Procedures (Enclose A)
SOURCE:	Mouse Ascites - Strain RH
FORM:	Frozen Liquid
PURIFICATION:	Ascites are sonicated and subjected to low speed centrifugation in buffer
INACTIVATION:	Each batch is inactivated by ultrasonication Note: Although inactivated, please handle with appropriate precautions.
APPLICATIONS:	Final working dilutions must be determined by end user.
STORAGE TEMPERATURE:	Frozen -70°C (avoid repeated thaws)
SHIPPING TEMPERATURE:	Frozen



ENCLOSE A

**TOXOPLASMA GONDII ANTIGEN
PRODUCTION FLOW SHEET**

Intraperitoneal inoculation of mice with suspension of Toxoplasma Gondii strain RH



Ascites growth



Ascites collection in bidistilled water



Sonication in physiological solution



Centrifugation/Supernatant collection

